

MyKETO®

PROTEIN MAIN DISHES



CHARACTERISTICS

- + High in protein-, high in fibre
- + Different flavours
- + Meal replacement contributes to weight loss
- + Meal replacement contributes to the maintenance of weight after weight loss
- + 25% recomanded daily intake of all vitamins and minerals in one portion
- + Contains vitamin D, Selenium and Zinc
- + Easy to prepare, just add water
- + Suitable for diabetics
- + Can be developed as a lactose free, gluten free or soya product.

Product specification

Category	Total diet or meal replacement for weight control
Recommended	For adults, teenagers
Country of origin	Czech Republic
Ingredients	Based on request for flavour and preferences. Pasta Bolognese: Ingredients: soy protein isolate, instant eggless noodles (wheat flour, vegetable fat, table salt, sugar, acidity regulators: sodium carbonate, potassium carbonate, thickener: guar gum), soy granules, tomato powder, soy flour, rice protein, inulin, potassium chloride, vegetable oils (sunflower, rapeseed), thickener: xanthan; dried beets, dried parsley, dried carrots, parsnips dried, dried onions, dried garlic, sodium chloride, antioxidant: rapeseed lecithin; acidity regulator: citric acid, spices, vitamin premix (maltodextrin, L-ascorbic acid, nicotinamide, DL- α -tocopheryl acetate, calcium D-pantothenate, riboflavin, pyridoxine hydrochloride, thiamine mononitrate, folic acid, D-biotin, cyanocobalamin), cholecalciferol, calcium phosphate, magnesium oxide, zinc lactate, manganese sulphate, copper sulphate, chromium picolinate, sodium iodide, sodium selenite, antioxidant: rosemary extract.

Instruction for use Pour the content of the sachet into a bowl, pour 140-160 ml of boiling water and mix well. Let swell for 4-5 minutes until the temperature of the soup is suitable for consumption.

	Packaging		Shelf life & Storage
Packaging	Sachet 1 serving (40-55 g)	Shelf life	12 months
	Metal/PE sachets	Storage conditions	Store in dry place (up to 60% relat. humidity) and at a temperature up to 25 °C.

BUSINESS ADDRESS: Michelska 18/12a | 140 00 Prague 4, Czech Republic
 FACTORY ADDRESS: V Lukach 923 | 392 01 Sobeslav, Czech Republic
 sales@goldim.cz | www.goldim.cz

